



SWIFTPET+

pipetting aid



always charged and ready for use
suitable for 0.5 ml to 100 ml pipettes
protected by filter and safety valve
ergonomically shaped handle
LCD display indicates the level
of battery charge



SWIFTPET+ is a revolutionary engine powered pipetting aid designed for cordless work with glass or plastic pipettes in the 0.5-100 ml range. Carefully modeled lightweight handle, together with smooth trigger buttons guarantee effortless pipetting even during extensive use. For convenient storage our device is equipped with "rest wings". Easily accessible switches allow to choose different operation modes depending on the volume of a pipette and viscosity of liquid. In order to protect the device against overfilling we equipped the **SWIFTPET+** with PTFE filters and safety valve, blocking any liquid from entering the unit. The filter and pipette holder can be easily exchanged and autoclaved. LCD display indicates setting options and level of battery charge.



Ordering information



Model	Cat. no.
SWIFTPET ⁺ with charging stand EU Plug	0351
SWIFTPET ⁺ with charging stand US Plug	0352
SWIFTPET ⁺ with charging stand UK Plug	0353
SWIFTPET ⁺ with charging stand Australia Plug	0354
Charging Stand	5410

Distributed by:

